



## ● Ion Torrent™ Technology

Run Validation and Introduction to the Torrent Browser  
(Torrent Suite™ Software v4.0.2)

● The world leader in serving science


# Agenda

---

- Best places to find information
  - Ion Community
  - Torrent Browser
- Run report Analysis & validation
  - Live demo
  - Slides as reference


# Ion Torrent™ Community

---



<http://ioncommunity.lifetechnologies.com/welcome>

# Ion Torrent™ Community



Ion Community

Search Ion Torrent Community

Applications

Datasets



Protocols

Products

Programs & Events

## See the clear picture with Ion RNA-Seq.

Reveal the diversity of the transcriptome.  
Discover what you've been missing >



### INTRO TO ion torrent

by *life* technologies™

Learn more about the Ion Torrent workflow, protocols, register for upcoming webinars and view recorded webinars.

Learn more >

Dear Ion Community,


We are happy to announce the roll-out of the all new Ion Community. The new site brings a number of great changes to the community.

- **New user interface:** The community look and feel has been redesigned from the ground up to bring a more modern look and feel to the site. The new primary navigation at the top allows for users to rapidly get to an area of interest in the community. The new design also leverages common elements to bring a more consistent look feel across all our software products.
- **Improved search functionality:** The process of identifying content of interest has been streamlined through visual layout (e.g., the primary navigation), organization of spaces, as well as the introduction of an improved search algorithm.

Ask For Help

Get Help

UPCOMING EVENTS

 07 **RNA-Seq Analysis Webinar**  
August 07 (Online)

# Ion Torrent™ Community

## INTRO TO ion torrent

by *life* technologies™

Learn more about the Ion Torrent workflow, protocols, register for upcoming webinars and view recorded webinars.

[Learn more >](#)

ion torrent

by *life* technologies™

Ion Community

Search Ion Torrent Community

Applications

Datasets

Protocols

Products

Programs & Events

Ion Torrent Community

Intro to Ion Torrent

Overview

Content


People

Subspaces and Projects

### INTRODUCTION TO ION TORRENT

Ion Torrent™ Semiconductor Sequencing  
Explore Ion Torrent™ Semiconductor Sequencing technology. Watch the video.

Ion Torrent™ next-gen sequencing technology



Register for Ion Torrent™ [Training](#) programs.

### RECENT ACTIVITY

Ion PGM™ Sequencer

Modified by Lifetech Administrator · 5 days ago

Expand

### SURVEY

What type of experiments are you working on?  
[Complete the survey and join the discussion.](#)

### RECENT CONTENT

Ion PGM™ Sequencer

Ion Proton Sequencer Throughput Competition

Intro to Ion Torrent

Ion Torrent - Everyone Deserves a Chance to Break Through

Transcriptome Webinar Series

More

### ACTIONS

Invite friends to join

Receive email notifications

Track in Communications

View feeds

# Ion Torrent™ Community

The screenshot displays the Ion Torrent Community website. At the top, the navigation bar includes the Ion Torrent logo, the text 'Ion Community', a search bar, and tabs for Applications, Datasets, Protocols, Products, and Programs & Events. The main content area is titled 'Ion Torrent Community → Products → Torrent Suite'. Below this, there are tabs for Overview, Content, People, Subspaces and Projects, and Reports. The 'Overview' tab is active, showing a 'WELCOME TO TORRENT SUITE SOFTWARE' message and a 'Power your analysis' section with a power button icon. A sidebar on the right lists 'ION SOFTWARE' products, with 'Torrent Suite™' circled in red. An arrow points from this circle to a 'GETTING STARTED' box, which contains links for 'Intro to Ion Torrent Informatics', 'Data Analysis Best Practices', 'Provide Product Feedback', and 'Training Programs'. Below this, a 'RECENT CONTENT' section lists various forum posts and documentation links.

ion torrent  
by life technologies™

Ion Community

Search Ion Torrent Community

Applications Datasets Protocols Products Programs & Events

ION SOFTWARE

ION PARTNERS

ION INSTRUMENTS

Ion Proton™ Sequencer

Ion PGM™ Sequencer

Ion AmpliSeq Designer™

Ion Reporter™

Torrent Browser Plugin Store™

Torrent Suite™

Product Feedback

Avadis NGS

BioTeam

DNASTar

Partek

SoftGenetics

Other Solutions

Overview Content People Subspaces and Projects Reports

Ion Torrent Community → Products → Torrent Suite

WELCOME TO TORRENT SUITE SOFTWARE

Power your analysis

Get the latest version of Torrent Suite™ Software

Learn about the latest Torrent Suite Software

Torrent Suite Documentation Links

RECENT ACTIVITY

TS analysis Error

nbahlis replied (in response to madhav)

Hi Madhav We did shut down our server and restart it with a monitor connected to it. On the reboot it stated that found 90 errors that I asked to fix (enter F). I restarted the analysis after this...

29 minutes ago

in Torrent Suite · 13 replies

Expand

GETTING STARTED

Torrent Suite

Intro to Ion Torrent Informatics

Data Analysis Best Practices

Provide Product Feedback

Training Programs

RECENT CONTENT

Filter by Categories & Tags

Re: TS analysis Error

Re: detailed sdk documentation

Torrent Variant Caller 3.6

VCF Tag Fields

Torrent Suite Documentation Links

More

# Ion Torrent™ Community

The screenshot displays the Ion Torrent Community website. At the top, the 'ion torrent by life technologies' logo is on the left, and a search bar labeled 'Search Ion Torrent Community' is on the right. Below the logo, a navigation bar contains links for 'Applications', 'Datasets', 'Protocols', 'Products', and 'Programs & Events'. The 'Protocols' link is highlighted, and a red circle is drawn around the 'ION PGM™ SEQUENCER' link in the dropdown menu. Below this, a green box contains links for 'Construct Library', 'Prepare Template', and 'Run Sequence' for both 'ION PGM™ SEQUENCER' and 'ION PROTON™ SEQUENCER'. The main content area shows the breadcrumb 'Ion Torrent Community → Protocols → PGM Protocols' and a sub-navigation bar with 'Overview', 'Content', 'People', 'Subspaces and Projects', and 'Reports'. The 'Overview' tab is selected, and a red circle is drawn around the 'ION PGM SYSTEM USER GUIDES AND PROTOCOLS' section. This section includes a question 'What kits should I use with the Ion PGM™ System? Review the decision trees:' followed by two links: 'Decision Tree for DNA Sequencing on the Ion PGM™ System' and 'Decision Tree for RNA Sequencing on the Ion PGM™ System'. Below these are three icons with corresponding links: a bar chart icon for 'Construct Library', a starburst icon for 'Prepare Template', and a DNA helix icon for 'Run Sequence'. To the right of this section is a 'GETTING STARTED' box with a 'Protocols' icon and links for 'Training Programs' and 'Provide Product Feedback'. Below that is an 'Ask For Help' section with a 'Get Help' button. At the bottom right is an 'ACTIONS' section with links for 'Start a discussion', 'Create by email', 'Invite friends to join', and 'Stop email notifications'. At the bottom left of the main content area, there are links for 'Torrent Suite Documentation' and 'PGM User Guide and Quick Reference Archive', and a link to 'Join the Ion PGM™ discussion'.

ion torrent  
by life technologies™

Ion Community

Search Ion Torrent Community

Applications Datasets Protocols Products Programs & Events

ION PGM™ SEQUENCER

ION PROTON™ SEQUENCER

Construct Library

Prepare Template

Run Sequence

Construct Library

Prepare Template

Run Sequence

Ion Torrent Community → Protocols → PGM Protocols

Overview Content People Subspaces and Projects Reports

ION PGM SYSTEM USER GUIDES AND PROTOCOLS

What kits should I use with the Ion PGM™ System? Review the decision trees:

- Decision Tree for DNA Sequencing on the Ion PGM™ System
- Decision Tree for RNA Sequencing on the Ion PGM™ System

Construct Library

Prepare Template

Run Sequence

GETTING STARTED

Protocols

Training Programs

Provide Product Feedback

Ask For Help

Get Help

ACTIONS

- Start a discussion
- Create by email
- Invite friends to join
- Stop email notifications

Torrent Suite Documentation

PGM User Guide and Quick Reference Archive

Join the Ion PGM™ discussion.



# Ion Torrent™ Community

The screenshot displays the Ion Torrent Community website interface. At the top, the 'Ion Community' header includes a search bar and navigation tabs for Applications, Datasets, Protocols, Products, and Programs & Events. The main content area is titled 'Ion Torrent Community → Products → PGM'. On the left, a sidebar contains links for Getting Started, User Guides and Bulletins, Tutorials and Training, FAQ, Run RecognitON Leader Board, and Ion Chip Recycling Program. The central 'GETTING STARTED' section includes links for PGM Training Programs and Product Feedback. The 'PRODUCT UPDATES' section features a list of Qubit® 2.0 Fluorometer Files, including various template kits, firmware, and calculators. A red circle highlights the 'Ion PGM™ Sequencer' link in the 'ION PARTNERS' section, with red arrows pointing to the 'GETTING STARTED' and 'PRODUCT UPDATES' sections.

**ION PARTNERS**

ION PARTNERS
Avadis NGS
BioTeam
DNAStar
Partek
SoftGenetics
Other Solutions

**GETTING STARTED**

PGM

- Training Programs
- Provide Product Feedback

**PRODUCT UPDATES**

**Qubit® 2.0 Fluorometer Files**

For information on using the following Qubit files, refer to the *Upgrade the Qubit® 2.0 Fluorometer firmware and software* section in either of the following User Guides:

- Ion PGM™ 200 Xpress™ Template Kit
- Ion OneTouch™ System
- Ion OneTouch 200 Template Kit v2
- Ion OneTouch™ 200 Template Kit v2 DL (for up to 200 base-read Sequencing)

Program (.bin) and upgrade (.qbt) files:

- Download: V3.10 Qubit\_FW\_MainCPU.bin
- Download: V3.10 Qubit\_FW\_UsbHost.bin
- Download: Ion\_PluginV310\_AF.qbt

USB Drives Tested for Use with Qubit® 2.0 Fluorometer


- Conversion Factor
- Qubit 2.0 Easy Calculator





# “Best Practices”, Ion Community


Check  
for  
updates!


## PRODUCT OVERVIEWS AND BEST PRACTICES (ION PGM SEQUENCER AND TORRENT SUITE)


 [Part I: Overview of Ion PGM™ System Protocols](#) →  
Guidelines for library and template preparation, sequencing,

 [Part II: Overview of Torrent Sequencing and Alignment \(v2\)](#)  
Sequence generation and alignment for Ion semiconductor d

 [Part III: Advanced Topics in Filtering of Ion Sequencing Data](#)  
Options to reanalyze Ion sequencing data for improved data

 [Part IV: Overview of Variant Calling in Ion Sequencing Data](#)  
Calling SNPs and indels in Ion semiconductor sequencing d


 [Part V: Advanced Topics in Variant Calling \(v3.4.x\)](#)  
Optimal sample-specific performance with the Torrent Varian

 **Part I: Overview of Ion PGM™ and Proton™ System Protocols**



 Version 24

created by [sue](#) on Aug 6, 2012 1:45 AM, last modified by [Shauna Clark](#) on Feb 15, 2013 2:18 PM

### Guidelines for Ion library and template preparation, sequencing, and instrument handling

This document outlines the best practices for successful Ion library construction, template preparation, and sequencing for Ion validated kits and protocols. The recommendations in this guide are only a brief summary of important protocol points, and first-time users should familiarize themselves with the appropriate user guides prior to beginning a protocol. User guides and quick references may be found in the  [Protocols](#) space on the [Ion Community](#).

### General guidance

We recommend using the latest library, template, and sequencing kits. For a list of the latest products and recommended kit compatibility, please see the  [Decision Tree for the Ion PGM™ System](#) or  [Decision Tree for the Ion Proton™ System](#). Users should also register for the [Ion Torrent Updates Group](#) on the Ion Community to receive email notifications regarding all Ion Torrent upgrades. For questions, contact technical support in your region or your local FAS.

### General laboratory guidance

- When designing the laboratory layout, consider the need for space separation of pre- and post-PCR activities. Separating the amplicon source, post-PCR activities from pre-PCR activities, and dedicating laboratory supplies and/or equipment to each space significantly reduces the potential for contamination

[http://ioncommunity.lifetechnologies.com/community/products/torrent\\_suite/best\\_practices](http://ioncommunity.lifetechnologies.com/community/products/torrent_suite/best_practices)

# TORRENT BROWSER™

Data analysis software to take you from primary to tertiary analysis, including optimized signal processing, base calling, sequence alignment, and variant analysis.



## PLAN

Plan instrument runs using Ion Torrent run templates or create your own.



## MONITOR

View the progress of the instrument and assess run metrics in real time.



## REVIEW

Search and review data across all runs, and drill down to see data in the run report.



## EXPORT

Automatically package and group data, and export it to other applications.

# New Account Registration - Request

ion torrent  
Sequencing for all™

Member Name  Password   ☐ Remember Me  
OR [Register](#) for a new account

1. Click Register
2. Fill out form
3. Submit
4. Approval

Iontraining

## New User Registration

Username:

E-mail:

Password:

Password (again):

New pending user registration for 'karenfriery'. Please visit [Account Management](#) (as an admin user) to review.

## Thanks for Registering

Your account has been created and is pending approval by an admin user.

Please contact your local admin and ask them to visit the Account Configuration page to activate your new account.

# New Account Registration – Admin Approval

**User Registrations** New user registrations awaiting approval

Username	Email	Full Name	Date Joined	
dleontest	darryl.leon@lifetech.com		Dec. 2, 2012, 2:46 p.m.	<input type="button" value="Approve"/> <input type="button" value="Reject"/>


1. Login as Admin
2. Review request
3. Approve/Reject

**Confirm Approve Registration** ×

Are you sure you want to approve the registration for **dleontest**?

\*Note: Only administrators can reset user passwords.

# Log into Torrent browser



## 3 main tabs

- **Plan:** Set up your run
- **Monitor:** Check your run in real time
- **Data:** Analyze your run report(s)

## Frederick Training Server

Plan

Monitor

Data

Plan Runs

Samples

Templates

Planned Run List

Search ...

New PGM Protocols/Manuals  
New Proton Protocols/Manuals  
IonCommunity  
Support  
**Local Documentation**  
System Support Diagnostics  
Instrument Diagnostics

## Torrent Suite Software Documentation

Welcome to the Torrent Suite Software documentation home page.

## Quick links to user documentation

- [Quickstart](#)
- [Torrent Browser User Interface Guide](#)
- [Torrent Browser Analysis Report Guide](#)
- [Use Cases](#)
- [Technical Notes and Whitepapers](#)

## Quick links to administrator documentation

- [Torrent Suite™ Software Data Management Guide](#)
- [Torrent Suite™ Software Administrator Documentation](#)

## Ion Community

For more help and the very latest information on Torrent Suite Software, check out Torrent Suite Software places on the Ion Community :


- [Torrent Suite Software](#)
- [Torrent Suite Software Best Practices](#)

## Icon

Access to the Ion Community requires (free) registration.

The Torrent Suite Software community is available to help, with resources and support, bringing together the power of shared ideas. In participating, you benefit from the collective knowledge and advice of a growing community of bioinformatics software developers.


After registering, use Torrent Suite Software resources to develop your own applications, furthering your own research projects and contributing to the community. There, you can view and download source code and samples, and receive notifications of the latest updates.

 [Torrent Suite™ Software User Documentation](#)

 [Torrent Suite™ Software Administrator Documentation](#)

Choose your  
documentation search  
method

# Planning a Run by Sample or Template





# Planning a Run by Sample

Frederick Training Server

Plan Monitor Data

Plan Runs Samples Templates Planned Run List

Plan Sample Type

Sample Sets

Import Samples from File Enter Samples Manually Sample Attributes

Sample Sets / Import Samples

Import Samples

1 : Select Sample File to Import :

Select file

Sample File (csv)

Sample File Format

Click button to download sample file format

2 : Select a Sample Set to Receive Samples :

Search sample set names Search Add Sample Set...

3 : Save & Finish

Uploads, parses & validates file, and saves if no errors found  
No samples imported yet

Summary

Create sample set(s) with samples in the imported sample sheet

Add Sample

Sample Name (required) : Required

FYI: Sample name in the sample set cannot be changed since this sample set has been planned.

Sample External ID : Optional

Gender : .....

Description : Optional

Role : .....


Relationship Group :

FYI: A unique number (e.g., 1) to designate multiples

Cancel Done

# Planning a Template Run

Create new run templates and plan instrument runs. You can use Ion Torrent templates or create new ones that match your own protocols.

 Whole-Genome Seq | [Add New Template](#) | [Plan New Run](#)

 RNA Seq | [Add New Template](#) | [Plan New Run](#)


 Ampliseq DNA | [Ampliseq.com Import](#) | [Add New Template](#) | [Plan New Run](#)


 Ampliseq Exome | [Ampliseq.com Import](#) | [Add New Template](#) | [Plan New Run](#)

 Ampliseq RNA | [Ampliseq.com Import](#) | [Add New Template](#) | [Plan New Run](#)

 TargetSeq | [Add New Template](#) | [Plan New Run](#)

 16S Target Sequencing | [Add New Template](#) | [Plan New Run](#)

 Generic Sequencing (Other) | [Add New Template](#) | [Plan New Run](#)



## Set up a run

- Use pre-installed template (application/product-based) or Create your own
- Start run on PGM™ sequencer using run code of planned run

# Review Settings of Pre-installed Template



Ampliseq DNA | [AmpliSeq.com Import](#) | [Add New Template](#) | [Plan New Run](#)

Ion AmpliSeq Cancer Panel 1_0 Lib Chem	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Inherited Disease Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Comprehensive Cancer Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Cancer Hotspot Panel v2	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>

Navigation: 1 2

### Review Template: Ion AmpliSeq Inherited Disease Panel

<b>Application</b> Application: AmpliSeq DNA	<b>Monitor</b> Bead Loading (%) : ≤ 30 Key Signal (1-100) : ≤ 30 Usable Sequence (%) : ≤ 30	<b>Plugins</b> Selected plugin(s): No selection saved with the Plan., Variant Frequency:
<b>Kits</b> Sample Preparation Kit: Library Kit Type: Ion AmpliSeq 2.0 Library Kit Forward Library Key: TCAG Forward 3' Adapter: ATCACCGACTGCCCATAGAGAGGCTGA GAC Templating Kit: Ion PGM Template OT2 200 Kit Sequencing Kit: IonPGM200Kit-v2	<b>Reference</b> Reference Library: hg19 Target Regions BED File: /hg19/unmerged/detail/4477686_IDP_b edfile_20120613.bed Hotspot Regions BED File:	<b>Projects</b> <b>Export</b> <b>Notes</b>
<b>Flows:</b> 500 Barcode Set (optional): Control Sequence (optional):		

Close

# Plan Run from a Template



[Ampliseq DNA](#) | [AmpliSeq.com Import](#) | [Add New Template](#) | [Plan New Run](#)

Ion AmpliSeq Cancer Panel 1_0 Lib Chem	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Inherited Disease Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Comprehensive Cancer Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Cancer Hotspot Panel v2	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>

6 - 9 of 9 items

Run Parameters

IonReporter

Application

Kits

Monitoring

Reference

Plugins

Output

Plan

Template Name: :

Ion AmpliSeq Inherited Disease Panel

Run Plan Name (required): :

Ion AmpliSeq Inherited Disease Panel

Enter a sample name for each plan (required at least one sample): :

Chip #	Sample Name (required)	Sample ID	Sample Description	Sample Tube Label
1	Your favorite sample			
2				

## Summary

Ion Reporter: None

Application: DNA

Sample Grouping:

Target Technique: AmpliSeq DNA

Workflow:

# Ion Reporter Integration

The screenshot displays the Ion Reporter web interface. At the top, there are tabs for 'Plan', 'Monitor', and 'Data'. Below these, a navigation bar includes 'Plan Runs', 'Samples', 'Templates', 'Planned Run List', and a highlighted button 'Create Plan from Ion AmpliSeq Inherited Disease Panel'. The main content area features a breadcrumb trail: 'Run Parameters' > 'IonReporter' > 'Application' > 'Kits' > 'Monitoring' > 'Reference' > 'Plugins' > 'Output' > 'Plan'. The 'IonReporter' step is highlighted with a red box. Below the breadcrumb, the text 'Select the IonReporter workflow to use:' is followed by three radio button options: 'None' (selected), 'CCPv1 from a Tumor/Normal Pair', and 'Genetic Disease Screening from an Exome Family Trio'. To the right of these options is a 'Sample Grouping' section with five radio button options: 'Sample\_Control', 'Self', 'Tumor\_Normal', 'Trio', and 'Other'. On the right side of the interface, a 'Summary' box displays the following information: 'Ion Reporter: None', 'Application: DNA', 'Sample Grouping:', 'Target Technique: AmpliSeq DNA', and 'Workflow:'. A blue 'Next ->' button is located at the bottom center of the main content area.

Plan Monitor Data

Plan Runs Samples Templates Planned Run List **Create Plan from Ion AmpliSeq Inherited Disease Panel**

Run Parameters **IonReporter** Application Kits Monitoring Reference Plugins Output Plan

Select the IonReporter workflow to use:

☒ None

☐ CCPv1 from a Tumor/Normal Pair

☐ Genetic Disease Screening from an Exome Family Trio

Sample Grouping

☐ Sample\_Control

☐ Self

☐ Tumor\_Normal

☐ Trio

☐ Other

Summary

Ion Reporter: None

Application: DNA

Sample Grouping:

Target Technique: AmpliSeq DNA

Workflow:

Next ->

Select the IR Workflow and Sample Grouping

# Template application

Plan

Monitor

Data

Plan Runs

Samples

Templates

Planned Run List

Create Plan from AmpliSeq Exome

Run Parameters

IonReporter

Application

Kits

Monitoring

Reference

Plugins

Output

Plan

Select the application, sample grouping and target technique, then hit next.

Application

Target Technique

☒ DNA

☐ RNA

☐ Metagenomics

☐ Typing

☐ Other

☐ AmpliSeq DNA

☐ TargetSeq

☐ Whole Genome

☒ AmpliSeq Exome

Previous

Next

Summary

Ion Reporter:

None

Application:

DNA

Sample Grouping:

Target Technique:

AmpliSeq Exome

Workflow:

# Plan Run from a Template - Kits

Run Parameters

IonReporter

Application

Kits

Monitoring

Reference

Plugins

Output

Plan

Select the sequencing kits and then hit next.

Sample Preparation Kit:

Library Kit Type: 

Details +

Ion AmpliSeq 2.0 Library Kit

Template Kit 

OneTouch

IonChef

 :  

Ion PGM Template OT2 200 Kit

Sequencing Kit:  

Ion PGM Sequencing 200 Kit v2

Flows :  

500

← Previous

Control Sequence (optional):

Chip Type (required):  

Ion 316™ Chip

Barcode Set (optional):

Mark as PCR Duplicates ☐ :

Next →

Summary

Ion Reporter: None

Application: DNA

Sample Grouping:

Target Technique: AmpliSeq DNA

Workflow:

Sample Preparation Kit:

Library Kit Type: Ion AmpliSeq 2.0 Library Kit

Library Key: Ion TCAG (TCAG)

Template Kit: Ion PGM Template OT2 200 Kit



# Plan Run from a Template - Monitor

Select the monitoring threshold preferences and then hit next.

Bead Loading (%): 30

Key Signal (1-100): 30

Usable Sequence (%): 30

← Previous

Next →

Summary

Ion Reporter: None

Application: DNA

Sample Grouping: AmpliSeq DNA

Target Technique: AmpliSeq DNA

Flow Transfer: 520 of 520 flows

(from Monitor tab)

73.6 % Loading

95.7 % Live ISPs

98.6 % Library ISPs

72 Key Signal

65 % Usable Seq

Set thresholds for run, view in Monitor Tab:

- Blue bar = greater than threshold set
- Red bar = less than threshold set

# Plan a Run from a Template - Reference

The screenshot shows a web interface for planning a run. At the top is a navigation bar with tabs: Run Parameters, IonReporter, Application, Kits, Monitoring, **Reference** (highlighted with a red box), Plugins, Output, and Plan. Below the navigation bar, the instruction reads: "Select the reference and regions to focus your analysis (optional), then hit next." There are three dropdown menus: "Reference Library" with "hg19(Homo sapiens hg19)" selected, "Target Regions File" (empty), and "Hotspot Regions File" (empty). At the bottom left is a "Previous" button, and at the bottom right is a "Next" button. On the right side, a "Summary" panel displays the following information:

Summary	
Ion Reporter:	None
Application:	DNA
Sample Grouping:	
Target Technique:	AmpliSeq DNA
Workflow:	
Sample Preparation Kit:	
Library Kit Type:	Ion AmpliSeq 2.0 Library Kit

**References**, sequences used in data analysis; Upload BED files

- Genome or other reference sequence
- Target Regions
- Hotspot Regions

# Plan a Run from a Template - Plugins

Select which plugins you want executed, then hit next.

Select All Clear Selections

<input type="checkbox"/> Assembler	<input type="checkbox"/> coverageAnalysis	<input type="checkbox"/> coverageAnalysisLite
<input type="checkbox"/> CreatFastqFile	<input type="checkbox"/> EnhancedQCwithPartekFlow	<input type="checkbox"/> ERCC_Analysis
<input type="checkbox"/> FastQC	<input type="checkbox"/> FastqCreator_3.6.2	<input checked="" type="checkbox"/> FileExporter
<input type="checkbox"/> PartekFlow	<input type="checkbox"/> TorrentSuiteCloud	<input checked="" type="checkbox"/> variantCaller <a href="#">Configure</a>
<input type="checkbox"/> variantCaller_3.6.2		

← Previous Next →

**Summary**

Ion Reporter: Brian 1.6

Application: DNA

Sample Grouping: AmpliSeq DNA

Target Technique: AmpliSeq DNA

Workflow: Sample Preparation Kit

**Plugins**, specific analysis invoked for this run

- Commercial plugins
- Customer plugins
- Plugins with parameter settings (somatic/germline)


# Plan Run from a Template - Output

The screenshot shows a software interface with a top navigation bar containing tabs: Run Parameters, IonReporter, Application, Kits, Monitoring, Reference, Plugins, Output (highlighted with a red box), and Plan. Below the tabs, a message reads: "Select the project(s) that will receive data from runs planned in this template, then hit next." A list of projects is shown with checkboxes: stuff, tm1, training (checked), TrainingClass, XLIN\_test, and YangTest. Below the list are a search bar, a "Search" button, and an "Add Project..." button. At the bottom left is a "← Previous" button, and at the bottom right is a "Next →" button. On the right side, a "Summary" panel displays the following information: Ion Reporter: None, Application: DNA, Sample Grouping: (blank), Target Technique: AmpliSeq DNA, Workflow: (blank), Sample Preparation Kit: (blank).

**Projects**, ways to mark data into groups

- Data sets can be in many projects
- Can group data from a cohort together
- Data in projects can be combined

# Plan Run from a Template



Run Parameters | IonReporter | Application | Kits | Monitoring | Reference | Plugins | Output | **Plan**

Template Name: :  
Ion AmpliSeq Inherited Disease Panel

Run Plan Name (required): :

Enter a sample name for each plan (required at least one sample): :

Chip #	Sample Name (required)	Sample ID	Sample Description	Sample Tube Label
1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

**Summary**

Ion Reporter: None

Application: DNA

Sample Grouping:

Target: AmpliSeq DNA

Technique:

Workflow:

## Plan Runs

- Accept template defaults and plan in seconds
- Add sample names to barcodes
- Plan for a group in one action

# Planned Runs List

Frederick Training Server

Plan Monitor Data

Plan Runs Samples Templates Planned Run List

### Planned Runs

All | by Template | by Sample

<input type="checkbox"/> Select	Run Code	Run Plan Name	Barcodes	Application	Project	Sample	Last Modified
<input type="checkbox"/>	53SJY	Ion AmpliSeq Inherited Disease Panel				None	2012/10/08 12:16 AM
		IonXpr					
		Panel_Yvonne2					
<input type="checkbox"/>	8EQW4	Ion AmpliSeq Comprehensive Cancer Panel_Yvonne					2012/10/05 06:58 AM
<input type="checkbox"/>	NB11M	Ion AmpliSeq Comprehensive Cancer Panel				Yvonne	2012/10/05 06:56 AM

Delete Selected Clear Selection


Type Run Code into sequencing instrument

or


Select Run Plan name on sequencing instrument

Review Edit Copy Delete

# Pulling in a Planned Run to the PGM™ Sequencer



Type in "Run Code" or  
Browse list of Planned  
Runs



✓ Verify planned run  
settings  
✓ Approve before starting  
run



# Easily Upload Multiple Plans - Overview

## Download Plan Multiple Runs File



[Ampliseq DNA](#) | [AmpliSeq.com Import](#) | [Add New Template](#) | [Plan New Run](#)

Ion AmpliSeq™ BRCA1 and BRCA2 Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">Delete</a>
IonXpressRunTemplate	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">Delete</a>
Ion AmpliSeq Custom ID	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Custom	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Cancer Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
<div><span>1</span> <span>2</span> <span>1 - 5 of 9 items</span></div>	

Review Uploaded Planned Runs

batchPlanning_2012_11_27_20_49_58.csv [Read-Only]																	
	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
1	Template	Plan name	Sample	Sample pr	Library kit	Templating	Control se	Sequence	Chip type	Flows	Bead load	Key signal	Usable se	Reference	Target regi	Hotspot re	Plugins
2	Ion AmpliSeq Fan_Nov212012				Ion Xpress	Ion OneTouch 200 Te	(200bp)	Ion PGM(tm)		500	30	30	30	hg19			
3	Ion AmpliSeq Fan_Nov212012				Ion Xpress	Ion OneTouch 200 Te	(200bp)	Ion PGM(tm)		500	30	30	30	hg19			
4	Ion AmpliSeq Fan_Nov212012				Ion Xpress	Ion OneTouch 200 Te	(200bp)	Ion PGM(tm)		500	30	30	30	hg19			
5																	

Modify File for Planning Multiple Runs

Select  
Plan  
Multiple



**Ampliseq DNA** | [AmpliSeq.com Import](#) | [Add New Template](#) | [Plan New Run](#)

Ion AmpliSeq Cancer Panel 1_0 Lib Chem	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Inherited Disease Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Comprehensive Cancer Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Cancer Hotspot Panel v2	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
<div><span>⏪</span> <span>⏩</span> <span>1</span> <span>2</span> <span>⏭</span> <span>⏮</span> <span>⏯</span></div> <div>6 - 9 of 9 items</div>	

Enter  
Number of  
Plans to  
Create

**Plan Runs from Template: Ion AmpliSeq Cancer Hotspot Panel v2** ×

**Number of plans to create(required)**


enter a whole number between 2 and 100

Number of plans to create offline

Cancel

Download CSV for batch planning

# Save CSV File and Edit



batchPlanning_2012_11_27_20_49_58.csv [Read-Only]																	
	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
1	Template	Plan name	Sample	Sample pr	Library kit	Templating	Control se	Sequence	Chip type	Flows	Bead load	Key signal	Usable se	Reference	Target regi	Hotspot re	Plugins
2	Ion AmpliSeq	Fan_Nov212012			Ion Xpress	Ion OneTouch 200	Te	(200bp)	Ion PGM(tm)	500	30	30	30	hg19			
3	Ion AmpliSeq	Fan_Nov212012			Ion Xpress	Ion OneTouch 200	Te	(200bp)	Ion PGM(tm)	500	30	30	30	hg19			
4	Ion AmpliSeq	Fan_Nov212012			Ion Xpress	Ion OneTouch 200	Te	(200bp)	Ion PGM(tm)	500	30	30	30	hg19			
5																	

# Upload Plans

## Templates



Favorites



Favorites

Upload Plans

TiN\_Round4

Review | Plan Run | Plan Multiple | Copy | Edit | Delete

Generic

Review | Plan Run | Plan Multiple | Copy | Edit | Delete

1

## Upload Plan Runs

CSV File for Plans

Browse...



Select CSV file to add

One row for each plan.

Tips:

(1) Plan values not included in the CSV file will default to values of the template you are planning from.

(2) Plugins and Project name, allow more than 1 attribute. Use ; as separator (e.g., enter My Project1; My Project2 for Project name)

Cancel

Upload CSV for batch planning

# Review uploaded Planned Runs

Plan

Monitor

Data







Plan Runs

Samples

Templates

Planned Run List

Planned Runs

<input type="checkbox"/> Select	Run Code	Run Plan Name ▲	Barcodes	Application	Project	Sample	Last Modified ▼	
<input type="checkbox"/>	6FI2A	A3				B3	2012/11/27 09:15 PM	
<input type="checkbox"/>	7KA5A	A2				B2	2012/11/27 09:15 PM	
<input type="checkbox"/>	74KTS	A1				B1	2012/11/27 09:15 PM	

# Plan cloning: Copy

The screenshot shows a web application interface with a top navigation bar containing 'Plan', 'Monitor', and 'Data' tabs. Below the tabs are links for 'Templates' and 'Planned Runs'. The main content area is titled 'Planned Runs' and contains a table with the following columns: Select, Run Code, Run Plan Name, Barcodes, Application, Project, Sample, and Last Modified. The table lists four planned runs. A context menu is open over the third row, which has 'Testing1256' as the Run Plan Name. The menu options are Review, Edit, Copy, and Delete. A red arrow points to the 'Copy' option.

Select	Run Code	Run Plan Name	Barcodes	Application	Project	Sample	Last Modified
<input type="checkbox"/>	PVZWK	darryl				good	2012/12/18 11:25 AM
<input type="checkbox"/>	4TKVU	Testing123	none			E20641-pool40-L1756	2012/12/18 11:25 AM
<input type="checkbox"/>	ML1MV	Testing1256	none			E20641-pool40-L1756	2012/12/18 11:25 AM
<input type="checkbox"/>	NWQFU	Testing1256	none			E20641-pool40-L1756	2012/11/15 07:11 AM


# Easily Create Run Template from Ion AmpliSeq™ Designer downloaded files

**Ion AmpliSeq™  
Designer**

The screenshot shows the 'Darryl Cancer Hotspot Panel' design results. It includes a table of gene targets and their coverage.

Type	Name	Chrom	Start	End	Target (bp)	Missed (bp)	Covered (%)
Gene (CDS Only)	PTEN1				1,212	17	98.6
Gene (CDS Only)	TP53				1,426	48	96.63
Gene (CDS Only)	LAS1				2,412	71	97.06

Additional information: IAD32385 - Results ready, 97.31% Coverage, 2 (20 ng) Pools (Input DNA), 200 bp Amplicon Size, 5.05 kb Target Size.



**Torrent  
Browser**

The screenshot shows the 'Ampliseq DNA' section of the Torrent Browser. It lists two designs: 'Darryl Cancer Hotspot Panel' and 'Ion AmpliSeq Custom ID'. Each design has links for 'Review', 'Plan Run', 'Plan Multiple', 'Copy', 'Edit', and 'Delete'.

Design Name	Actions
Darryl Cancer Hotspot Panel	Review   Plan Run   Plan Multiple   Copy   Edit   Delete
Ion AmpliSeq Custom ID	Review   Plan Run   Plan Multiple   Copy



# Download zip file from Ion AmpliSeq™ Designer to local computer

The screenshot shows the Ion Torrent by Life Technologies web interface. The user is logged in as 'Hi, Darryl Leon' and is viewing the 'Darryl Cancer Hotspot Panel'. The design 'mynewdesign1' is selected. The results for 'IAD32385' are ready. The design parameters are: 97.31% Coverage, 2 (20 ng) Pools (Input DNA), 200 bp Amplicon Size, and 5.05 kb Target Size. The 'Download results' button is highlighted with a red box. A red arrow points from this button to a file explorer window showing the file '4475346\_CHP2\_resultsD.zip' on the desktop.

ion torrent by life technologies™

Hi, Darryl Leon | Sign Out

Dashboard My Designs Panels My Account Help

Darryl Cancer Hotspot Panel

Switch design: mynewdesign1 Edit Copy Add design

IAD32385 - Results ready More solutions +

Ion AmpliSeq™ Library Kit 2.0  
Your design requires 2 library kit reactions per sample

97.31% Coverage ? 2 (20 ng) Pools (Input DNA) ? 200 bp Amplicon Size ? 5.05 kb Target Size

Order design Download results More actions +

End	Target (bp)	Missed (bp)	Covered (%)	
	1,212	17	98.6	UCSC
	1,426	48	96.63	UCSC
	2,412	71	97.06	UCSC

< Previous Showing 1-3 of 3 results Next > Last >> 15 Per page

Open


Look in: Desktop

My Recent Documents Desktop My Documents My Computer My Network Places

V2.2\_TorrentSuite V3.0\_TorrentSuite V3.2.1\_TorrentSuite V3.2\_TorrentSuite V3.4\_TorrentSuite V3.6\_Torrent Suite V3.8\_Torrent Suite WebExInfo WeeklyStatus 4475346\_CHP2\_resultsD.zip AmpliSeqOrderinDemoandTSRNaanalysis.txt basecaller\_filter.xlsx

File name: Files of type: All Files Open as read-only Open Cancel

# Navigate to Ion AmpliSeq™ Designer Zipped file and Upload



# Verify Zipped BED File is Added

---

---

## Available BED Files

Upload BED files

## Upload a BED file for hg19:


4475346\_CHP2\_resultsD.zip (225 KB) 100%

Select a new BED File





Hotspot ☐ ?

Upload

# Verify Template is Automatically Generated


 **Ampliseq DNA** | [AmpliSeq.com Import](#) | [Add New Template](#) | [Plan New Run](#)

Darryl Cancer Hotspot Panel	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">Delete</a>
Ion AmpliSeq Custom ID	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Custom	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>
Ion AmpliSeq Cancer Panel Template Example	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">Delete</a>
barcoded template	<a href="#">Review</a>   <a href="#">Plan Run</a>   <a href="#">Plan Multiple</a>   <a href="#">Copy</a>   <a href="#">Edit</a>   <a href="#">Delete</a>

  **1** [2](#) [3](#)  

1 - 5 of 14 items

# Select References in Torrent Browser



# Select Reference Genome

PlanMonitorData

AboutReferencesServicesPluginsConfigureAccounts

Reference Sequences

Obsolete Reference Sequences

Target Regions

Hotspots

Test Fragments

Barcodes

Upload History

Reference Sequences

Get Ion ReferencesAdd Reference Sequence

Name	Description	Notes	Date ▼	Index version	Status	Enabled
hg19	Homo sapiens hg19	AmpliSeq supported hg19	2013/09/18 03:01 PM	tmap-f3	found	true
PRRS	Porcine Respiratory and Reproductive Syndrome Virus VR2332		2013/08/29 11:08 PM	tmap-f3	created	true
jathropaCP	JCchloroplast	jatropha chloroplast	2013/08/29 05:01 PM	tmap-f3	created	true
jatcp	jatropha	Jatropha curcas chloroplast, complete genome	2013/08/29 03:48 PM	tmap-f3	created	true
Linnocua	Listeria innocua		2013/08/29 03:47 PM	tmap-f3	created	true
mitoJSS	mitoRef	test mito	2013/07/26 11:11 AM	tmap-f3	created	true
dm_cdna_fa_zip	dm cdna fa zip	none	2013/06/26 12:06 PM	tmap-f3	created	true
C_diff	C diff	OpGen Reference	2013/03/28 08:41 PM	tmap-f3	created	true
rod	Chrodney		2013/03/12 02:12 PM	tmap-f3	created	true
am	chrn1_AM		2013/03/12 02:12 PM	tmap-f3	created	true
hg19_rna	Homo sapiens RNA Canonical Transcript	AmpliSeq supported hg19	2012/08/06 04:07 PM	tmap-f3	downloading	true
mm10	Mus musculus 10	AmpliSeq supported mm10	2012/08/06 04:07 PM	tmap-f3	created	true

# Target File Selection

PlanMonitorData

AboutReferencesServicesPluginsConfigureAccounts

## New Target Regions

Target Regions File :

Select File

Please select a BED file to upload.

Reference :

A\_hydrophila\_ATCC7966 - Aeromon

Description :

optional

Notes :

optional

Upload Target Regions File

Cancel

# Monitor Tab

Frederick Training Lab Hi, ionadmin Help Sign Out

**Plan** **Monitor** **Data**

Runs in Progress

**Runs in Progress** Live Update On

List View | [Table View](#)

**user\_ROB-74**  
Robin

[3.2\\_Re-analysis](#)

**Base Calling** No Run Plan

Run Started: 2012/07/17 12:42 PM  
Chip: 314  
Notes:

Flow Transfer: 520 of 520 flows

77.4 %  
Loading

91.3 %  
Live ISPs

98.6 %  
Library ISPs

?  
Signal

? %  
Usable Seq

Information is populated in real time



# Data Tab – Table View & List View

Frederick Training Lab

Hi, ionadmin Help Sign Out

Plan Monitor **Data**

Completed Runs & Results Projects

Completed Runs & Results

List View **Table View** Live Update On

★ □ Date

Any Sample

Any Instrument

**Table View** is useful for scanning, searching, or sorting a large number of results

Run Name	Sample	Appli...	Run Date	Analysis Date	Status	Chip	Rep Report Name	Ref Genom	Barcode	Flov	T...	M...	Q...	O...
☆ user_ROB-74	am2		2012/...	2012/...	Base ...	314	<a href="#">3.2_Re-analysis</a>	e_c...						
☆ user_BAT-7...	BAT_FRDCI...		2012/...	2012/...	Com...	314	<a href="#">3.2_Re-analysis</a>	e_c...		47...	127	54...	60...	
☆ user_CAT-11	Jenny-lindsay		2012/...	2012/...	Com...	314	<a href="#">3.2_Re-analysis</a>	e_c...		65...	236	13...	15...	
☆ user_ROB-75	jun912		2012/...	2012/...	Com...	314	<a href="#">3.2_Reanalysis</a>	e_c...		45...	215	80...	97...	
☆ user_ROB-76	ryan913		2012/...	2012/...	Com...	314	<a href="#">ROB-76_ReAna_TS3.0</a>	e_c...		57...	240	11...	13...	
☆ user_BAT-83	Ecoli_1to13...		2012/...	2012/...	Com...	314	<a href="#">Auto_BAT-83_170</a>	e_c...		34...	230	64...	79...	

Frederick Training Lab

Hi, ionadmin Help Sign Out

Plan Monitor **Data**

Completed Runs & Results Projects

Completed Runs & Results


List View **Table View** Live Update On

★ □ Date

Any Sample

Any Instrument


**List View** provides more detail at a glance; useful for reviewing a run and all of its reports at once

☆  **user\_ROB-74**  
Robin  
Flows Complete 2012/07/17 12:42 PM

Reference: e\_coli\_dh10b  
Sample: am2  
Chip: 314  
Project: [Class](#)

Reanalyze Edit Delete

Report Name	Total Reads	Mean Read Len.	Q20 Bases	Output	Date	Status
☆ <a href="#">3.2_Re-analysis</a>					2012/10/07 11:56 PM	Alignment
☆ <a href="#">Auto_ROB-74_159</a>	571 k	237	112 M	135 M	2012/08/16 01:18 PM	Completed

☆  **user\_BAT-77-BAT\_FRDClass\_1\_130**  
Batman

Reference: e\_coli\_dh10b  
Sample: BAT\_FRDClass\_1\_130  
Chip: 314

Reanalyze Edit

# Data Tab

Frederick Training Lab

Hi, ionadmin Help Sign Out

Plan Monitor **Data**

Completed Runs & Results Projects

**Search** by run name, run date, chip type, etc.

Completed Runs & Results

List View **Table View**

★ ☐ Date Search names Go Any project

Sample Any Reference All Flows Any Chip

Instrument All Result Status Results new to old Clear

**Table of Runs & Run Reports**

	Name	Sample	Appli...	Run Date	Analysis Date	Status	Chip	Rep Report Name	Ref Genom	Barcode	Flov	T... R...	M... R... L...	Q... B...	O...
☆	user_ROB-74	am2		2012/...	2012/...	Base ...	314	<a href="#">3.2_Re-analysis</a>	e_c...	...	...	...	...	...	...
☆	user_BAT-7...	BAT_FRDCI...		2012/...	2012/...	Com...	314	<a href="#">3.2_Re-analysis</a>	e_c...	...	47...	127	54...	60...	...
☆	user_CAT-11	jenny-lindsay		2012/...	2012/...	Com...	314	<a href="#">3.2_Re-analysis</a>	e_c...	...	65...	236	13...	15...	...
☆	user_ROB-75	jun912		2012/...	2012/...	Com...	314	<a href="#">3.2_Reanalysis</a>	e_c...	...	45...	215	80...	97...	...
☆	user_ROB-76	ryan913		2012/...	2012/...	Com...	314	<a href="#">ROB_76_ReAna_TS3.0</a>	e_c...	...	57...	240	11...	13...	...
☆	user_BAT-83	Ecoli_1to13...		2012/...	2012/...	Com...	314	<a href="#">_BAT-83_170</a>	e_c...	...	34...	230	64...	79...	...

Click to see run report


Reanalyze  
Edit  
Add Report to Project

Open a report

# Run Report - Overview



- Navigation bar {
- Barcoding table {  
(if applicable)
- Unaligned/  
Aligned metrics {
- Plugin summaries {
- Output files {



# Run Report – Top

The screenshot displays the 'Frederick Training Lab' web interface. At the top right, there is a user status bar with 'Hi, ionuser', a 'Help' dropdown, and a 'Sign Out' link. Below this is a navigation bar with three main tabs: 'Plan', 'Monitor', and 'Data'. The 'Data' tab is highlighted with a red box. Below the 'Data' tab, there is a sub-navigation bar with 'Completed Runs & Results' (highlighted with a red box) and 'Projects'. Under 'Completed Runs & Results', there are three sub-tabs: 'Run Summary' (highlighted with a red box), 'Output Files', and 'Plugin Summary'. A red callout box labeled 'Navigate run report' points to the 'Run Summary' tab. Below the sub-navigation bar, the main content area shows 'Run Summary: R\_2012\_09\_13\_12\_29\_55\_user\_CAT-11'. To the right of this, there is a 'Reports' dropdown menu. The dropdown is open, showing three options: '3.2 Re-analysis (241)' (highlighted with a red box), '3.0 Re-analysis (239)', and 'Auto\_CAT-11\_171 (236)'. A red callout box labeled 'See other run reports for this run\*' points to the dropdown menu.

Frederick Training Lab

Hi, ionuser Help Sign Out

Plan Monitor Data

Completed Runs & Results Projects

Run Summary Output Files Plugin Summary

Reanalyze Summary PDF Plugins PDF Classic Report

Run Summary: R\_2012\_09\_13\_12\_29\_55\_user\_CAT-11






Reports 3.2 Re-analysis (241) 3.0 Re-analysis (239) Auto\_CAT-11\_171 (236)

Navigate run report

See other run reports for this run\*

- \* Most recent analysis defaults as representative run report. User can change representative run report.


## Run Report – Barcoding information (if used)

Barcode Name	Sample	Bases	>=Q20 Bases	Reads	Mean Read Length	Read Length Histogram	Files
No barcode	None	1,855,688	1,516,254	25,940	71 bp		<a href="#">BAM</a> <a href="#">BAI</a>
IonXpress_001	None	29,580	24,580	268	110 bp		<a href="#">BAM</a> <a href="#">BAI</a>
IonXpress_002	None	24,070	20,293	219	109 bp		<a href="#">BAM</a> <a href="#">BAI</a>
IonXpress_003	None	54,162	44,683	485	111 bp		<a href="#">BAM</a> <a href="#">BAI</a>
IonXpress_004	None	1,927,550	1,580,007	14,833	129 bp		<a href="#">BAM</a> <a href="#">BAI</a>

For each Barcode, we show...

- Total output (bases)
- AQ20 Bases
- AQ20 reads
- Mean AQ20 read length

# Run Report – Summary of Unaligned reads




Efficiency of chip loading

Percentage of wells that generated reads that passed filtering

Average read lengths


Information about bases and reads before mapping

# Chip Loading Image: Very Good Loading Example





Ion 314™v2 chip heat map, indicating percent loading across the physical surface

ISP Density






# Run Report validation – Key Signal



The **Consensus Key 1-Mer** graph shows the strength of the signal from the first three one-mer bases of the library key. This graph represents the consensus signal measurement of release of H<sup>+</sup> during nucleotide incorporation. The y-axis shows signal strength, measured in **Counts**, which is an arbitrary but consistent unit of measure. The x-axis shows time as nucleotide **Flows** over the chip.



# Run Report – Chip Well Details




- Measure loading performance and enrichment efficiency
- Quantify read filter metrics: assess polyclonality and low quality reads

# Ion Sphere™ Particle (ISP) Identification Summary

<b>Addressable Wells</b>	<b>1,262,519</b>	
With ISPs	1,100,543	87.2%
Live	1,020,656	92.7%
Test Fragment	10,853	01.1%
Library	1,009,803	98.9%
<b>Library ISPs</b>	<b>1,009,803</b>	
Filtered: Polyclonal	296,525	29.4%
Filtered: Low Quality	53,864	05.3%
Filtered: Primer Dimer	2	00.0%
<b>Final Library ISPs</b>	<b>659,412</b>	<b>65.3%</b>

Percentage of sequence available for analysis after filtering

Mouse over for help



# Ion Sphere™ Particle (ISP) Identification Enrichment Calculation

<b>Addressable Wells</b>	<b>1,262,519</b>	
With ISPs	1,100,543	87.2%
Live	1,020,656	92.7%
Test Fragment	10,853	01.1%
Library	1,009,803	98.9%
<b>Library ISPs</b>	<b>1,009,803</b>	
Filtered: Polyclonal	296,525	29.4%
Filtered: Low Quality	53,864	05.3%
Filtered: Primer Dimer	2	00.0%
<b>Final Library ISPs</b>	<b>659,412</b>	<b>65.3%</b>

% enrich =

Lib ISPs

1,009,803


( Wells with ISPs

1,100,543

- TF ISPs


10,853

# Run Report – Read Length Details



Histogram of read length for all reads in output files

Read Length Details



# What does Q20 mean?

---

## Predicted Quality

Q: What does Q17/Q20 mean?

A: The **Q20 length** value attributed to a read is an estimate of the read length at which the predicted total error rate in the read will correspond to a Phred-scale quality score of 20. The Phred scale is defined as  **$-10 \times \log_{10}(\text{error probability})$** , so Q20 corresponds to an error rate of 1% and Q17 corresponds to an error rate of 2%. Please, consult the following link for a more complete description of Phred values: [http://en.wikipedia.org/wiki/Phred\\_quality\\_score](http://en.wikipedia.org/wiki/Phred_quality_score).

The **Q20 length** is determined by looking at the per-base quality scores to estimate the total read error rate at every position in the read. For example, if the first base had a predicted quality score of 20 and the second base had a predicted quality score of 17, then the estimated total read error for the first two bases would be  **$0.5 \times 1\% + 0.5 \times 2\% = 1.5\%$** . This total read error rate is evaluated at every position in the read, after which the maximal position at which the total read error rate is 1% or less is identified. This position defines the Q20 length of the read.

The Q20 length is derived entirely from the predicted per-base quality scores. The predicted quality scores are somewhat conservative and in cases where alignment to a reference sequence is also possible users will generally find that the predicted Q17 lengths tend on average to be shorter than the corresponding AQ17 lengths which are based on the actual as opposed to predicted errors (see the following section on AQ20/AQ17)


# What does AQ20 mean?

---

- Alignment of reads is a process to map the read to a reference
- Any discrepancy in alignment to a reference is listed as a mismatch, and described by a mapping quality score
- In Torrent Browser, throughput reads are  $\geq$  AQ20 reads
- AQ20 length of a read is the longest length at which the aligned read has a mismatch rate of 1% or less
- AQ20 content is the sum of all the AQ20 lengths

# Run Report – Reads Aligned to Reference

Aligned to E. coli DH10B



Shows metrics gathered after mapping to reference.

# Run Report – Output Files

Hi, ionuser Help Sign Out

Reanalyze Summary PDF Plugins PDF Classic Report

Reports 3.2\_Re-analysis (241)

3.0\_Re-analysis | Frederick Training Lab - Torrent Browser Page 1 of 5

Run: R\_2012\_09\_13\_12\_29\_55\_user\_CAT-11

Run Summary

Unaligned

155M 71 660

86% Read Loading

ISP Density

3.0\_Re-analysis | Frederick Training Lab - Torrent Browser Page 2 of 5

Aligned to E. coli chr100

155M 33.2X 99.1%

Plugin Output

Test Fragments

Test: chr100 chr100 chr100 chr100 chr100 chr100 chr100 chr100 chr100 chr100

Read length histogram

100 200 300 400 500 600 700 800 900 1000

File format of older version

## Output Files

File Type	Unaligned Reads	Aligned Reads
Library	BAM SFF FASTQ	BAM BAI
Barcodes	BAM SFF FASTQ	BAM BAI

Barcode Name	Sample	Bases	>=Q20 Bases	Reads	Mean Read Length	Read Length Histogram	Files
No barcode	None	3,919,092	3,098,417	42,495	92 bp		BAM BAI
IonXpress_007	None	809,623,089	671,266,731	4,673,571	173 bp		BAM BAI

10 items per page 1 - 2 of 2 items

File formats produced for extensive data analysis via third party software packages (offered on Webstore) or other bioinformatics tools.



# Plugin Summary



## Select a plugin

- [Alignment](#) — v3.6.56201
- [AmpliconCoveragePlots](#) — v1.5.1
- [Assembler](#) — v3.4.2.0
- [combineAlignments](#) — v2.2.3-31037
- [coverageAnalysis](#) — v3.6.58977
- [ERCC\\_Analysis](#) — v3.6.56201
- [FastQC](#) — v3.4.1.1
- [FastqCreator](#) — v3.6.0-r57238
- [FilterDuplicates](#) — v3.6.58660
- [IonReporterUploader](#) — v3.6.0-r58488
- [IonReporterUploader\\_V1\\_2](#) — v3.2.0-r50766
- [RunRecognition](#) — v3.6.56201
- [sampleID](#) — v3.6.56201
- [SAMStat](#) — v3.0.0
- [SFFCreator](#) — v3.6.56201
- [TorrentSuiteCloud](#) — v3.6.58782
- [variantCaller](#) — v3.6.59049
- [variantCaller\\_3.4.2](#) — v3.4.51874

## variantCaller Features

Supports TargetSeq™, Ion AmpliSeq™, & whole genome applications

- Each may be run as Germ Line or Somatic

Enables calling for all custom products

- Supports any genome or custom BED file

Summarizes statistics for barcodes, coverage, chromosomes, variants & hotspots


Calls at 5% or 50%

Compatible with barcoded samples

Planning page for automatic variant calling

# Test Fragment Summary

Plugin Summary **Test Fragments** Analysis Details Support Software Version H+

Test Fragment	Reads	Percent 50AQ17	Read length histogram
TF_A	56,898	98%	

⏪ ⏴ 1 ⏵ ⏩ 1 - 1 of 1 items

# Run Report Validation: Test fragment Summary

TF\_A

```
TGTTTTAGGGTCCCCGGGGTTAAAAAGGTTTCGAACTCAACAGCTGTCTGG  
CAGCTCGCTCTACGATCTGAGACTGCCAAGGCACACAAGGGATAGG
```

- Test Fragment designed with variable homopolymeric challenges (2mer, 3mer and 4mer)
- TFs are used for install validation
- TFs are spiked into sample for run validation
- No specification on test fragment sequencing performance
- Specific kit type may have different TFs

# Report Information

<div> <div>Plugin Summary</div> <div>Test Fragments</div> <div>Analysis Details</div> <div>Support</div> <div>Software Version</div> </div>	
Run Name	R_2013_11_15_09_36_24_user_BAT-120
Run Date	Nov. 15, 2013, 9:36 a.m.
Run Flows	500
Projects	<a href="#">training</a>
Sample	e.coli
Sample Tube Label	None
Reference	e_coli_dh10b
Instrument	Batman
Flow Order	TACGTACGTCTGAGCATCGATCG
Library Key	TCAG
TF Key	ATCG
Chip Check	Passed
Chip Type	314R
Chip Data	single
Barcode Set	
Analysis Name	Auto_user_BAT-120_456
Analysis Date	Nov. 15, 2013, 2:38 p.m.
Analysis Flows	500
runID	LXTUK

<div> <div>Plugin Summary</div> <div>Test Fragments</div> <div>Analysis Details</div> <div>Support</div> <div>Software Version</div> </div>	
Torrent_Suite	4.0.1
host	68JKKQ1
ion-analysis	4.0.5-1
ion-dbreports	4.0.21-1
ion-gpu	4.0.0-1
ion-pipeline	4.0.6-1
ion-plugins	4.0.21-1
ion-torrentr	4.0.4-1
Script	21.5.5
LiveView	545
DataCollect	462
OS	20
Graphics	35

# Customer Support Archive (CSA)

## Output Files

File Type	Reads	Aligned Reads
Library	<span>BAM</span> <span>SFF</span> <span>FASTQ</span>	<span>BAM</span> <span>BAI</span>

Plugin Summary Test Fragments Analysis Details **Support** Software Version

Select plugins to run + Expand All - Collapse All Refresh plugin status

- Download the **Customer Support Archive**
- [View the report log](#)

Files for Troubleshooting;  
useful for Technical Support

If you need troubleshooting /support, send in your CSA file (have it ready!)

# About Link

**Frederick Training Server**

Plan Monitor Data

About References Services Plugins Configure Accounts

**Releases**

Product	Version
Torrent Suite	4.0.2
• ion-analysis	• 4.0.6-1
• ion-dbreports	• 4.0.39-1

About  
References  
Services  
Plugins  
Configure  
Accounts

## More Information and Assistance

- [IonCommunity](#)
- [Support](#)
- [Local Documentation](#)
- [System Support Diagnostics](#)
- [Instrument Diagnostics](#)

# Services Link

## Frederick Training Server

[Plan](#)[Monitor](#)[Data](#)

[About](#)[References](#)[Services](#)[Plugins](#)[Configure](#)[Accounts](#)

Jobs Server

Hostname  
68JJKQ1

Service Name  
RSM\_Launch  
RabbitMQ  
celery\_periodic  
celery\_plugins

Configuration

Enabled	Auto-Acknowledge	Signal Processing Input	Basecalling Input	Output Files	Intermediate Files
<a href="#">Configure</a> <input type="checkbox"/>	Disabled	Disabled	Disabled	Disabled	Archive

Data Archive

/media/external

Mounted

Mounted

[Configuration Log](#)[History](#)

### Disk Usage

Filesystem	Used Space	Used Space	Size (GB)	Free (GB)	Raw Data (GB)	Reports (GB)	Other (GB)
File Server Space							
/results/	<div><div>Used</div><div>free</div></div> <div>59.4%</div>	10837.2	4395.8	5151.4	605.0	134.5	
Archive: /media/external							
/media/external	<div><div>Used</div><div>free</div></div> <div>76.5%</div>	7.5	1.8				

# Configure Link

Frederick Training Lab

Hi, ionadmin Help Sign Out

Plan Monitor Data

About References Services Plugins **Configure** Accounts

About  
References  
Services  
Plugins  
**Configure**  
Accounts

## Customer Support Contact

This is the person in your organization who should be notified during a support request of problems related to the nature of an experiment/run.

## IT Contact

This is the person in your organization who should be notified during a support request of problems related to the Torrent Server hardware or the network environment.

### Customize Site Name

Frederick Training Lab [Save Name](#)

### Email

[Add Email](#)

Edit	Enabled	Email Address
<a href="#">Edit</a>   <a href="#">Delete</a>	<input type="checkbox"/>	LT8@gmail.com
<a href="#">Edit</a>   <a href="#">Delete</a>	<input type="checkbox"/>	LT7@gmail.com
<a href="#">Edit</a>   <a href="#">Delete</a>	<input type="checkbox"/>	LT6@gmail.com
<a href="#">Edit</a>   <a href="#">Delete</a>	<input checked="" type="checkbox"/>	LT5@gmail.com
<a href="#">Edit</a>   <a href="#">Delete</a>	<input checked="" type="checkbox"/>	LT3@gmail.com

### Database Administration

The [Admin Interface](#) provides direct access to the database entries for system administrators.



# Accounts Link

Frederick Training Lab

Hi, ionadminHelpSign Out

PlanMonitorData

AboutReferencesServicesPluginsConfigureAccounts

Settings

AboutReferencesServicesPluginsConfigureAccounts

**User Profile / Account Information** Allows to view and update your existing user profile and account information

Usernameionadmin

Full Name[Not Set]

REST API Key

eb2b0c33fb6d0f3c1fc2c4951adad80e7cbfad46

Account LevelAdmin

Name

Email

ionadmin@iontorrent.com

Phone number

ResetSubmit

**User Registrations** New user registrations awaiting approval

No user accounts pending approval.

---

# ion torrent



by *life* technologies™

For Research Use Only. Not for use in diagnostic procedures.

Life Technologies is a Thermo Fisher Scientific brand. © 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries.