

Best Practices: Workflow QC

The world leader in serving science

#### **Best Practices**

Focus on General Ion Workflow and Importance of QC Checkpoints:

- Starting material
  - Quality and Quantity
- Library Construction
  - Post-Shear Profile
  - Final Library Profile
  - Quantification
- Templated Bead Preparation
  - Qubit® → *Pre*-enrichment
- Sequencing
  - Run Reports



#### "Best Practices", Ion Community



## PRODUCT OVERVIEWS AND BEST PRACTICES (ION PGM SEQUENCER AND TORRENT SUITE)

 Part I: Overview of Ion PGM™ and Proton™ System Protocols

ii Version 24

created by sue on Aug 6, 2012 1:45 AM, last modified by Shauna Clark on Feb 15, 2013 2:18 PM

- Part II: Overview of Torrent Sequencing and Alignment (v2.: Sequence generation and alignment for Ion semiconductor da
- Part III: Advanced Topics in Filtering of Ion Sequencing Day Options to reanalyze Ion sequencing data for improved data of
- Part IV: Overview of Variant Calling in Ion Sequencing Data Calling SNPs and indels in Ion semiconductor sequencing da
- Part V: Advanced Topics in Variant Calling (v3.4.x)
  Optimal sample-specific performance with the Torrent Variant

## Guidelines for Ion library and template preparation, sequencing, and instrument handling

This document outlines the best practices for successful Ion library construction, template preparation, and sequencing for Ion validated kits and protocols. The recommendations in this guide are only a brief summary of important protocol points, and first-time users should familiarize themselves with the appropriate user guides prior to beginning a protocol. User guides and quick references may be found in the Protocols space on the Ion Community.

#### General guidance

We recommend using the latest library, template, and sequencing kits. For a list of the latest products and recommended kit compatibility, please see the ☐ Decision Tree for the lon PGM™ System or ☐ Decision Tree for the lon Proton™ System. Users should also register for the lon ☐ Torrent Updates Group on the lon Community to receive email notifications regarding all lon Torrent upgrades. For questions, contact technical support in your region or your local FAS.


#### General laboratory guidance

 When designing the laboratory layout, consider the need for space separation of pre- and post-PCR activities. Separating the amplicon source, post-PCR activities from pre-PCR activities, and dedicating laboratory supplies and/or equipment to each space significantly reduces the potential for contamination

http://ioncommunity.lifetechnologies.com/community/products/torrent\_suite/best\_practices



#### Ion Workflow Overview



#### How QC Fits in Workflow


Reminder: All examples herein are for gDNA as starting material unless otherwise noted

**Templated** Sequence + Starting Library Bead Run Construction Material = DNA Evaluation Preparation **□**Quantity □Calculate TDF □Post-Shearing □Check Run **Profile** Report **□**Quality □Qubit® Assay: □Check size post-amplification, selection unenriched sample □Quantify library

Last but not least....

Do not mix kit components and align with the most current protocol

#### A Closer Examination



#### Starting Material: Qualify

- Ion Xpress<sup>™</sup> Plus gDNA Fragment Library Preparation
- Begin with pure and intact starting material


#### Appendix A. Evaluate the quality of the genomic DNA

#### Assess the integrity and size by gel electrophoresis

We recommend checking the integrity and size of your DNA preparation by gel electrophoresis. Use of a spectrophotometer to assess DNA quality can be misleading, because many molecules absorb in ultraviolet light.


#### Figure 2 Examples of genomic DNA preparations

Examples of high-quality DNA with no contaminating RNA (lanes 1 and 2), compared to lower quality samples containing RNA contamination (lanes 3 and 4). The RNA runs as a diffuse smear at the bottom of the gel. M is a lambda *HindIII* molecular weight marker.



### Starting Material: NanoDrop® Considerations

#### **Example: optimal NanoDrop® curve**




- ■260/230 ratio = Organic chemicals and solvent contamination Important to have this ratio above 1.6 for QPCR
- ■260/280 ratio = Protein Contamination Ratios between 1.8 – 2.0 are considered pure

DNA has no ethanol contamination, curve is normal shape with OD260 having higher absorbance than OD230




#### Starting Material: NanoDrop® Considerations

#### **Example: Residual Ethanol**




- ■OD 230 reading is much higher than OD 260
- ■This may result in inaccurate concentration readings and affect downstream enzymatic steps, especially when DNA concentration is low

#### Qubit® 2.0 Fluorometer vs. NanoDrop® Spectrophotometer



Qubit<sup>®</sup> Instrument: effective range 10pg/μL to 1μg/μL DNA NanoDrop<sup>®</sup> Instrument: effective range 2ng/μL to 15μg/μL


#### A Closer Examination




#### Post Shearing Profile

#### Shearing: Example of Ion Shear™ Plus Reagent digestion

- Applications (starting material) may drive profile of sheared DNA
- Shearing target size is read length dependent
- Shearing conditions may need to be optimized (e.g., incubation time and temp. with enzymatic methods)
- Input amount → expected recovery after size selection





#### Final Library Profile

#### Analyze size distribution

What are we looking for?

- Does observed size match expected size?
- How much non-specific product exists relative to desired library?
  - Little or no remaining primer, primer-dimer, artifacts, by-products
- Options
  - Bioanalyzer® (strongly recommended) or equivalent instrument
    - Qualitative and quantitative measure of library
    - Define regions for proper sizing and quantitation
    - Important to follow maintenance recommendations
  - Gels
- Requires more material
- Gives a qualitative output of size distribution




### Final Library Profile: Bioanalyzer® Traces


For optimal sequencing performance, libraries need to be appropriately sized before template preparation

Case

- 100 bp template prep = ~130 bp insert, ~200 bp library
- 200 bp template prep = ~260 bp insert, ~330 bp library

#### Fragment DNA libraries







Study

### Final Library Profile: Bioanalyzer® Traces


#### Whole transcriptome library (RNA)



#### smRNA library (RNA)



#### Ion AmpliSeq™ DNA library (Cancer Primer Pool)



## Final Library: Quantification prior to Template Preparation


- Quantification critical to success of template prep
  - To properly determine the molar concentration and library yield, consider the following:
    - Precision and accuracy of the measurement system
    - Sensitivity and specificity of the detection method
    - Stay within the quantitative range of the assay
- Measurement systems
  - Real-time PCR
  - Bioanalyzer® analysis (also gives library size distribution)




- Qubit® fluorometer
- Nanodrop® Instrument- Not recommended
  - "total nucleic acid" is measured: Detects anything at A260 (DNA, RNA, primers, dNTPs, etc)



## Final Library: Quantitative Range



- Quantitative range varies by assay
  - Note RFUs for marker, ladder, and sample peaks



Overloaded



#### A Closer Examination

Starting Material

Library Construction Templated
Bead
Preparation

Sequence +
Run
Evaluation

- ✓ Qualify
- ✓ Quantify
- ✓ Post-Shearing Profile
- ✓ Final Library Profile
- ✓ Quantify Library

- □Calculate TDF
- □Qubit® Assay: post-amplification (unenriched) sample

#### Templated Bead Prep QC

- To maximize system throughput:
  - Properly qualify and quantify your library
  - Calculate appropriate template dilution factor for your library type and template prep kit
  - Run the Qubit® assay (Ion Sphere Quality Control Kit) after the amplification, prior to enrichment.
    - Save 2 uL of your sample and store at 4°C if you skip this step.
       You can choose to run the assay later, if needed.
    - Use the correct Conversion Factor for each kit lot, available on the lon Community
  - Following guidance in the Template Prep User Guide, proceed to enrichment OR re-do step(s) as needed

#### Templated Bead Prep QC: Steps 1 & 2

- √1. Properly qualify and quantify your library
  - 2. Calculate appropriate template dilution factor for your library type and template prep kit

#### **Library Dilution Guidelines:**


Starting material	Input type	1x TDF	
DNA	gDNA - whole genome	~26 pM	
	amplicon – long		
	amplicon – short		
	Ion AmpliSeq v2.0	16-20 pM	
	fusion PCR	~26 pM	
RNA	WT: polyA-selected RNA	~32 pM	
	WT: rRNA-		
	small RNA	~32 pM	




Note: Depending on version of Ion OneTouch™ kit, volume of dilution added to aqueous phase varies

#### Templated Bead Prep QC: Step 3

3. Run the Qubit® assay (Ion Sphere™ Quality Control Kit ) after the template amplification, prior to enrichment



 The Percent Templated ISPs calculates automatically and is displayed for each sample:



#### Templated Bead Prep QC: Step 4

The recommended optimal range is not intended to be a pass/fail criteria. The range provides guidance for the quality of the sample.

Note: If the results are outside the desired Percent Templated ISPs range, then increase or decrease the library input appropriately.

Percent Templated ISPs	Description  Sample contains an insufficient number of templated ISPs to achieve optimal loading density on the Ion Chip.	
<10%		
10-30%	Optimal amount of library.	
>30%	Sample will yield multi-templated ISPs (mixed reads).	

### Templated Bead Prep QC Considerations

Sample	Qubit® 2.0 Fluorometer Phenotype	Ion Personal Genome Machine™ (PGM™) Possible Observation	Possible Root Cause	Recommended Action
	<10% Templated ISPs	Lower loading     Lower % enriched     Lower key signal     Lower throughput	Too little library input into template preparation	Increase library input to target 20–25% templated ISPs or     Continue with sequencing; expect lower throughput
	>30% Templated ISPs but <70%	Increased number of filtered reads	Too much library input into template preparation	Decrease library input to target 20–25% templated ISPs or     Continue with sequencing; expect lower throughput
Unenriched	>70% Templated ISPs	Low loading     Low % enriched     Lower throughput     High % filtered reads	Ion OneTouch™ Instrument underperformance	Troubleshoot with Life Technologies Technical Support or a Field Application Scientist
	primer filtered	Increased % primer dimer filtered reads     Lower throughput	Adapter dimer contaminating library, more likely in short amplicon, AmpliSeq <sup>™</sup> or miRNA libraries	Check bioanalyzer traces for adapter dimer peak [Amplicon tibrary or Ion AmpliSeq™ Library peak around 70 bp; miRNA Library peak around 60bp].     Re-purify Library using Agencourt® AMPure® XP Kit clean-up steps as outlined in the appropriate user guides.

#### A Closer Examination

## Starting Material

Library Construction Templated
Bead
Preparation

Sequence +
Run
Evaluation


- ✓ Qualify
- ✓ Quantify
- ✓ Post-Shearing Profile
- ✓ Final Library Profile
- ✓ Quantify Library

- √ Calculate TDF
- ✓ Qubit® Assay: post-amplification (unenriched) sample

□Check Run Report




#### Check Run Report




- ☐ Library Summary
- □ Test Fragments
- ☐ ISP Summary
- □ Software Version

The Run Report is a .pdf document that is *included* in the Customer Support Archive (CSA) file


#### Run Report: Library



#### Run Report: Test Fragments



#### Run Report: ISP Summary




#### Run Report: Software Version

ftware Version	Is the software version current?
Torrent_Suite	2.2
Datacollect	210
Graphics	18
UveView	345
os	19
Script	18.1.6
host	tahiti19
ion-alignment	2.2.4-1
ion-analysis	2.2.12-1
ion-gpu	1.2-1
ion-pipeline	2.2.12-1
Ion-tomentR	2.2.8-1

- Different software versions produce different "expected" outputs
- Some software versions are not compatible with certain applications...make sure you are using the most recent version

#### Bringing it all together



Starting
Material = DNA

Library Construction Templated
Bead
Preparation

Sequencing + Run Reports

- ✓ Post-Shearing Profile
- ✓ Final Library Profile
- Quantify Library


STOP!

Qubit®, qPCR, Bioanalyzer®:

Always inspect results

Example: Bioanalyzer® → ladder and markers are integral to accurate sizing and quantification

## Bioanalyzer® Instrument: Importance of Ladder & Markers



**After Correction**: [library] = **8.2 nM** (14 fold difference)

Starting Material


Library Construction Templated
Bead
Preparation

Run Report

	Count	Percentage
Library ISPs / Percent Enrichment	162843	52%
• Filtered: Too short	678	<1%
• Filtered: Keypass failure	7905	5%
Filtered: Mixed / Polyclonal	149905	92%
► Filtered: Poor Signal Profile	4264	3%
→ Final Library Reads	91	<1%

Note: Screen captures are from Torrent Suite™ v.1.4 software. Later versions of Torrent Suite™ software may look different





- Very low throughput: only 62 reads aligned
- Key signal is very good, both intensity and shape

Note: Screen captures are from Torrent Suite™ v.1.4 software. Later versions of Torrent Suite™ software may look different



#### Ion Sphere™ Particle (ISP) Identification Summary Percentage Count Loading Density (Avg ~ 26%) 100 % Total Addressable Wells 1262519 90 % Wells with ISPs 333985 26% 80 % Height = 1152 wells ---> 70 % Live ISPs 182178 55% 60 % Test Fragment ISPs 19335 11% 50 % Library ISPs 162843 89% 40 % 30 % Count Percentage 20 % Library ISPs / Percent Enrichment 162843 52% 10 % • Filtered: Too short 678 <1% <--- Width = 1280 wells ---> · Filtered: Keypass failure 7905 5% Filtered: Mixed / Polyclonal 149905 92% 4264 3% Filtered: Poor Signal Profile Final Library Reads 91 <1%

- Low loading density
- Extremely high % of reads rejected by mixed filter → Possible reasons?

Note: Screen captures are from Torrent Suite™ v.1.4 software. Later versions of Torrent Suite™ software look different



Starting Material

Library Construction Templated
Bead
Preparation

Run Report



Calculate TDF

Be cognizant of your dilution scheme (Follow Template Prep User Guide for guidance):



- Typically, high polyclonality indicates too much library was added into Templated Bead Preparation
- In this example, undiluted stock library was added into the amplification mix

- "Polyclonality happens"
  - Qubit % template-positive ISPs indirectly measures polyclonality
  - Higher %positive on Qubit means higher % polyclonal
  - Can change polyclonality by re-doing template prep with different amount of library

	Count	Percentage
Library ISPs / Percent Enrichment	630,569	84%
Filtered: Polyclonal	173,663	28%
Filtered: Primer dimer	75	<1%
Filtered: Low quality	90,938	14%
Final Library Reads	365,893	58%

	Count	Percentage
Library ISPs / Percent Enrichment	987,953	95%
Filtered: Polyclonal	228,637	23%
• Filtered: Primer dimer	223	<1%
► Filtered: Low quality	77,111	8%
→ Final Library Reads	681,982	69%

Consider % final library reads



#### Remember...

- QC Checkpoints are there for a reason
- Each is integral to optimizing PGM<sup>TM</sup> performance

# 

Life Technologies is a Thermo Fisher Scientific brand. © 2014 Thermo Fisher Scientific, Inc. All rights reserved. NanoDrop is a registered trademark of NanoDrop Technologies, LLC. Bioanalyzer is a registered trademark of Agilent Technologies, Inc. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.